

Development of Standards for Calibrated Flow Cytometry

Jacob Beal, NIST SBSC Flow Cytometry WG

COMBINE
Salt Lake City
Oct. 2015


NIST **Raytheon**
BBN Technologies

Overview

- **Need:** absolute unit measurements from large numbers of single cells
- **Technology:** Calibrated Flow Cytometry
- **Process:** NIST SBSC Flow Cytometry WG

Why single cell measurements?

Are these groups of cells the same?




But they've got the same mean fluorescence...?

Why large numbers of samples?


Example: RNA Replicon Cotransfection

Scatter of Red vs. Blue for ../042613/042613/042613_042613_26.fcs



Example: Constitutive mKate in HEK293

Color Compensation Model




Quantifying variation components requires per-cell measurements of large populations of cells

Why absolute units?

- Metrology 101: precision measurement enables:
 - Comparison of results across experiments and labs
 - Deeper insight into the behavior of devices
 - Effective dissemination of materials and methods
 - Testing and validation of materials and systems
 - Establishment of commercial & industrial services
 - Safety assurance and traceability of responsibility

Many biological measurements are only relative!

How Flow Cytometry Works




Challenges:

- Autofluorescence
- Variation in measurements
- Spectral overlap
- Time Contamination
- Lots of data points!
- Different protein fluorescence
- Individual cells behave (very) differently

Metrology vs. Flow Cytometry

Flow cytometry great for per-cell measurements, but...




Arbitrary unit output depends on...

- Instrument brand, configuration
- Interference from other colors
- Choice of instrument settings
- Run-to-run calibration drift

Fortunately, this can be corrected...

Calibrated Flow Cytometry



*Negative control:
autofluorescence*


Arbitrary Units (a.u.)

Net a.u.


Compensated a.u.

FITC a.u.


Calibrated MEFL



*Single positive controls:
compensate for spectral overlap*




*Multi-color control:
conversion of non-FITC channels*




*Calibration beads:
Convert FITC a.u. to MEFL*

Example: Predicting Repressor Cascades


Precision dose-response measurement allows high-precision prediction with quantitative models




Prediction of Repressor Cascade



Range vs. Error for 6 Cascades



How much does calibration matter?



NIST Synthetic Biology Standards Consortium

- Flow cytometry working group:
 - Initial goal: two standards documents
 - Minimum Standard for Flow Cytometry Calibration Information *what needs to be done*
 - Methods for Calibrated Flow Cytometry *a way to do it*
 - Process & Timeline:
 - Currently: drafting in progress
 - Nov 3: NIST SBSC meeting in San Francisco
 - Dec '15: Dissemination for critique & improvement
 - Feb '16: Version 1 ready for publication process
 - Onward: calibrate FSC, SSC, interlab studies

Join this open group: <https://groups.google.com/forum/#!forum/sb-sc-flow-cytometry>